

BLOOD BANK REFRIGERATORS AND FREEZERS FOR BLOOD AND PLASMA STORAGE

EMOBANK

MODELS	EMOBANK										
	100/2	200/3	200/4	300/5	300/6	700/5	700/6	700/7	700/8		
DIMENSIONS											
Internal dim. (WxDxH)	mm	480x500x450	480x500x815	480x500x1385				590x675x1500			
External dim. (WxDxH)	mm	600x615x850	600x615x1425	600x615x1995				740x815x1998			
Liters	l	95	176	303				600			
Weight	kg	85	100	115				170			
TECHNICAL SPECIFICATIONS											
Temperature range	°C	0 - +15									
Power voltage	V	230									
Noise level	dB	<50									
Max. absorbed power	A	1.5	2						1		
Heat dissipation	kcal/h	170	215						180		
Refrigerant gas		R 290						R290			
STRUCTURE											
External finish		Prepainted galvanized-plasticized steel sheet - color gray RAL 9006									
Internal material		AISI 304 Stainless Steel									
Insulation	mm										
Shelves/Drawers	N.	2	3	4	5	6	5	6	7	8	
Maximum load per shelf	kg	30					40				
Capacity per drawer/shelf	N.	30 bags (250ml), 16 bags (450ml)						76 sacche (250ml), 44 sacche (450ml)			
Maximum load	N.	60 bags (250ml), 32 bags (450ml)	90 bags (250ml), 48 bags (450ml)	120 bags (250ml), 64 bags (450ml)	150 bags (250ml), 80 bags (450ml)	180 bags (250ml), 96 bags (450ml)	380 bags (250ml), 220 bags (450ml)	458 bags (250ml), 264 bags (450ml)	532 bags (250ml), 308 bags (450ml)	608 bags (250ml), 352 bags (450ml)	

MODELS	EMOBANK						
	700/B	900/B	900/9	900/B	1500/12	1500/14	1500/B
DIMENSIONS							
Internal dim. (WxDxH)	mm	590x675x1500	660x875x1510			1330x675x1500	
External dim. (WxDxH)	mm	740x815x1998	810x1010x1998			1480x815x1998	
Liters	l	600	870			1400	
Weight	kg	170	175			230	
TECHNICAL SPECIFICATIONS							
Temperature range	°C	0 - +15					
Power voltage	V	230					
Noise level	dB	<50	<50			<50	
Max. absorbed power	A	1	1,3			1,2	
Heat dissipation	kcal/h	180	140			220	
Refrigerant gas		R290					
STRUCTURE							
External finish		Prepainted galvanized-plasticized steel sheet - color gray RAL 9006					
Internal material		AISI 304 Stainless Steel					
Insulation	mm	75					
Shelves/Drawers	N.	7	8	9	12	14	
Maximum load per shelf	kg	40					
Capacity per drawer/shelf	N.	76 bags (250ml), 44 bags (450ml)	90 bags (250ml), 55 bags (450ml)			76 bags (250ml), 44 bags (450ml)	
Maximum load	N.	532 bags (250ml), 308 bags (450ml)	720 bags (250ml), 440 bags (450ml)	810 bags (250ml), 495 bags (450ml)	912 bags (250ml), 528 bags (450ml)	1064 bags (250ml), 616 bags (450ml)	