

BLOOD BANK REFRIGERATORS AND FREEZERS FOR BLOOD AND PLASMA STORAGE

EMOBANK

MODELS	EMOBANK										
	100/2	200/3	200/4	300/5	300/6	700/5	700/6	700/7	700/8		
DIMENSIONS											
Internal dim. (WxDxH)	mm	480x500x450	480x500x815			480x500x1385		590x675x1500			
External dim. (WxDxH)	mm	600x815x850	600x815x1425			600x815x1995		740x815x1998			
Liters	l	108	196			332		584			
Weight	kg	85	100			115		170			
TECHNICAL SPECIFICATIONS											
Temperature range	°C	0 - +15									
Power voltage	V	230									
Noise level	dB	55	55			55		38,5			
Max. absorbed power	A	1,5	2			1					
Heat dissipation	heat/h	170	215			180					
Refrigerant gas		R 290			R 290		R290				
STRUCTURE											
External finish		Prepainted galvanized-plasticized steel sheet - color gray RAL 9006									
Internal material		AISI 304 Stainless Steel									
Insulation	mm										
Shelves/Drawers	N.	2	3	4	5	6	5	6	7	8	
Maximum load per shelf	kg	30					40				
Capacity per drawer/shelf	N.	30 bags (250ml), 16 bags (450ml)					76 sacche (250ml), 44 sacche (450ml)				
Maximum load	N.	60 bags (250ml), 32 bags (450ml)	90 bags (250ml), 48 bags (450ml)	120 bags (250ml), 64 bags (450ml)	150 bags (250ml), 80 bags (450ml)	180 bags (250ml), 96 bags (450ml)	380 bags (250ml), 220 bags (450ml)	456 bags (250ml), 264 bags (450ml)	532 bags (250ml), 308 bags (450ml)	608 bags (250ml), 352 bags (450ml)	
MODELS	EMOBANK										
	700/B	900/B	900/B	1500/12	1500/14	1500/B					
DIMENSIONS											
Internal dim. (WxDxH)	mm	590x675x1500	680x875x1510		1330x675x1500						
External dim. (WxDxH)	mm	740x815x1998	810x1010x1998		1480x815x1998						
Liters	l	584	812		1347						
Weight	kg	170	175		230						
TECHNICAL SPECIFICATIONS											
Temperature range	°C	0 - +15									
Power voltage	V	230									
Noise level	dB	38,5	43,8		43,8						
Max. absorbed power	A	1	1,3		1,2						
Heat dissipation	heat/h	180	140		220						
Refrigerant gas		R290									
STRUCTURE											
External finish		Prepainted galvanized-plasticized steel sheet - color gray RAL 9006									
Internal material		AISI 304 Stainless Steel									
Insulation	mm	75									
Shelves/Drawers	N.	7	8	12		14					
Maximum load per shelf	kg	40									
Capacity per drawer/shelf	N.	76 bags (250ml), 44 bags (450ml)	90 bags (250ml), 55 bags (450ml)		76 bags (250ml), 44 bags (450ml)						
Maximum load	N.	532 bags (250ml), 308 bags (450ml)	720 bags (250ml), 440 bags (450ml)	912 bags (250ml), 528 bags (450ml)		1084 bags (250ml), 616 bags (450ml)					