

STABILITY TEST AND PLANT GROWTH CHAMBERS

EKOCT - EKOCTL - EKOCH - EKOCHL

Model	EKOCH				EKOCT	
	700	1500	900P	1500P	700	1500
Comercial code	12880	12881	13379	13383	12871	12873
Capacity (l)	700	1500	900	1500	700	1500
Temperature range (°C)	-10°C / +60°C		-20°C / +60°C		+10°C / +50°C	
Dewpoint (°C)	+5°C				/	
Relative humidity (%)	20/90				/	
Humidity precision (%)	-/+4				/	
Supply voltage (V)	230v(+6 / -10%) 50-60 Hz 1+T(G)					
Maximum absorbed current (A)	11				8	9
Cooling gas	R134A		R452A		R134A	
External dimensions (mm)	740(L)x940(P)x2020(H)	1480(L)x940(P)x2020(H)	810(L)x1135(P)x2020(H)	1480(L)x940(P)x2020(H)	740(L)x815(P)x2020(H)	1480(L)x815(P)x2020(H)
Internal dimensions (mm)	590(L)x600(P)x1500(H)	1330(L)x600(P)x1500(H)	660(L)x800(P)x1500(H)	1330(L)x600(P)x1500(H)	590(L)x600(P)x1500(H)	1330(L)x600(P)x1500(H)
Overall weight with standard packaging (Kg)	180	270	210	270	180	270
Noise level (dbA)	< 65					

Model	EKOCHL			
	700	1500	900P	1500P
Comercial code	12884	12885	13413	13302
Capacity (l)	700	1500	900	1500
Temperature range (°C) with lights on	+10°C / +50°C			
Temperature range (°C) with lights off	+10°C / +50°C		-10°C / +60°C	
Dewpoint (°C)	+5°C			
Relative humidity (%)	20/90			
Intensity light (lux)	0 / 20.000			
Humidity precision (%)	+ 5			
Supply voltage (V)	230v(+6 / -10%) 50-60 Hz 1+T(G)			
Maximum absorbed current (A)	20			
Cooling gas	R134A		R452A	
External dimensions (mm)	840(L)x940(P)x2100(H)	1580(L)x940(P)x2100(H)	910(L)x1135(P)x2100(H)	1580(L)x940(P)x2100(H)
Internal dimensions (mm)	590(L)x600(P)x1500(H)	1330(L)x600(P)x1500(H)	660(L)x800(P)x1500(H)	1330(L)x600(P)x1500(H)
Overall weight with standard packaging (Kg)	190	280	220	280
Noise level (dbA)	< 65			

Model	EKOCTL			
	700	1500	900P	1500P
Comercial code	12882	12883	13414	13300
Capacity (l)	700	1500	900	1500
Temperature range (°C) with lights on	+10°C / +50°C			
Temperature range (°C) with lights off	+10°C / +50°C		-10°C / +60°C	
Intensity light (lux)	0 / 20.000			
Supply voltage (V)	230v(+6 / -10%) 50-60 Hz 1+T(G)			
Maximum absorbed current (A)	18			
Cooling gas	R134A		R452A	
External dimensions (mm)	840(L)x815(P)x2100(H)	1580(L)x815(P)x2100(H)	910(L)x1010(P)x2100(H)	1580(L)x815(P)x2100(H)
Internal dimensions (mm)	590(L)x600(P)x1500(H)	1330(L)x600(P)x1500(H)	660(L)x800(P)x1500(H)	1330(L)x600(P)x1500(H)
Overall weight with standard packaging (Kg)	190	280	220	280
Noise level (dbA)	< 65			