

BLOOD BANK REFRIGERATORS AND FREEZERS FOR BLOOD AND PLASMA STORAGE

EMOBANK

MODELS	EMOBANK									
	100/2	200/3	200/4	300/5	300/6	700/5	700/6	700/7	700/8	
DIMENSIONS										
Internal dim. (LxDxH)	mm	480x500x450	480x500x815		480x500x1385		590x 675x 1500			
External dim. (LxDxH)	mm	600x615x850	600x500x815	600x615x1425	600x615x1995		740x 815x 1998			
Liters	l	95	176		303		600			
Weight	kg	480x500x450	100		115		170			
TECHNICAL SPECIFICATIONS										
Temperature range	°C	0 - +15								
Climatic class		I								
Power voltage	V	230								
Noise level	dB	<50								
Max. absorbed power	A	1,5	2				4			
Heat dissipation	kcal/h	170	215				180			
Refrigerant gas		R 290		R 290				R290		
STRUCTURE										
External finish		Prepainted galvanized-plasticized steel sheet - color gray RAL 9006								
Internal material		AISI 304 Stainless Steel								
Insulation	mm									
Shelves	Nr.	2	3	4	5	6	5	6	7	8
Maximum load per shelf	kg	30				40				
Capacity per drawer/shelf	Nr.	30 bags (250ml), 16 bags (450ml)				76 bags (250ml), 44 bags (450ml)				
Maximum load	Nr.	60 sacche (250ml), 32 bags (450ml)	90 bags (250ml), 48 bags (450ml)	120 bags (250ml), 64 bags (450ml)	150 bags (250ml), 80 bags (450ml)	180 bags (250ml), 96 bags (450ml)	380 bags (250ml), 220 bags (450ml)	456 bags (250ml), 264 bags (450ml)	532 bags (250ml), 308 bags (450ml)	608 bags (250ml), 352 bags (450ml)

MODELS	EMOBANK						
	700/8	900/8	900/9	900/B	1500/12	1500/14	1500/B
DIMENSIONS							
Internal dim. (LxDxH)	mm	590x675x1500	660x875x1510		1330x675x1500		
External dim. (LxDxH)	mm	740x815x1998	810x1010x1998		1480x815x1998		
Liters	l	600	870		1400		
Weight	kg	170	180		230		
TECHNICAL SPECIFICATIONS							
Temperature range	°C	0 - +15					
Climatic class		I					
Power voltage	V	230					
Noise level	dB	<50	<60		<50		
Max. absorbed power	A	4			5		
Heat dissipation	kcal/h	180	140		260		
Refrigerant gas		R290					
STRUCTURE							
External finish		Prepainted galvanized-plasticized steel sheet - color gray RAL 9006					
Internal material		AISI 304 Stainless Steel					
Insulation	mm	75					
Shelves	Nr.	7	8	9	12	14	
Maximum load per shelf	kg	40					
Capacity per drawer/shelf	Nr.	76 bags (250ml), 44 bags (450ml)	90 bags (250ml), 55 bags (450ml)		76 bags (250ml), 44 bags (450ml)		
Maximum load	Nr.	532 bags (250ml), 308 bags (450ml)	720 bags (250ml), 440 bags (450ml)	810 bags (250ml), 486 bags (450ml)	912 bags (250ml), 528 bags (450ml)	1064 bags (250ml), 616 bags (450ml)	