

## PLASMA BLAST FREEZERS

# PLASMAFROST ITEM

Models	PlasmaFrost iTeM 3	PlasmaFrost iTeM 4
Internal shelves (n°)	3	4
Bags capacity (450ml)	24	40
Working temperature (°C)	-75	-75
Supply voltage (V)	230V (+6/-10%) 50 Hz/1+N+T(G)	400V (+6/-10%) 50 Hz/3+N+T(G)
Maximum absorbed current (A)	12	12
External dimensions (mm)(W,D,H)	1300x1490x1830	1300x1490x1830
Net weight (kg)	660	690
Noise level (dBA)	62	62